



Nucleus Hydro Dissection Cannula  
Chung-G.27  
B# 15-270



Nucleus Hydro Dissection Cannula  
G.25-26-27  
B# 15-255



Air Injection Cannula  
G.24-25-26-27-30  
B# 10-250



Double I / A Cannula (reverse)  
G.22 C# 20-115  
G.23 C# 20-135



Double I / A Cannula (original)  
G.22 B# 20-110  
G.23 B# 20-130



Iachrymal Irrigating Cannula  
B# 40-280



Mc Intyre Iachrymal probing Cannula  
22 ga x 50 mm Long Tube.Blunt tip  
23 ga x 50 mm Long Tube.Blunt tip  
With dual 0,3mm side opening  
B# 40-100  
B# 40-200



Silicon injection Cannula  
G16-18-19  
B#25-100



Double I & A Cannula  
J-Shaped  
B#20-116



Fluid Irrigating Cannula  
G.20  
F# 10-200



Sub Retinal Fluid Drainage  
Soft Tip G.20  
F# 10-95



Irrigating Cannula (Curved)  
For Lasik  
G.25  
B# 50-220 1 Hole  
B# 50-220/2 2 Hole  
B# 50-220/3 3 Hole  
B# 50-220/4 4 Hole  
B# 50-220/5 8 Hole



**Jensen**  
 Posterior Capsule Polisher  
 23 GA-Curved tube  
 Olive tip with Front opening port  
 B# 10-150




**Aspirating Cannula**  
 G.20  
 F# 10-85



**Simcoe**  
 cortex Extractor & Capsule Polisher  
 22&23 GA-Curved tube  
 sandblasted tip, port facing into Curve  
 B# 10-160



**Kratz**  
 Posterior Extractor & Capsule Polisher  
 23 GA-Curved tube  
 Sandblasted tip, port Facing into Curve  
 B# 10-170



**Self-Retaining Anterior Chamber Maintainer**  
 G.20(0.9mm) - B#30-300  
 G.21(0.8mm) - B#30-310  
 G.22(0.7mm) - B#30-320  
 G.23(0.6mm) - B#30-330  
 scale : 3/4



**Vitreotomy Infusion Cannula**  
 Infusion Cannula  
 6 mm  
 4 mm  
 B# 30-200  
 scale : 3/4



Bimanual Hand Pieces (Phaco)  
Integrated  
F# 10-151  
scale : 1



Bimanual Hand Pieces  
Hand Piece (phaco)  
F# 10-150  
scale : 1