



Bowman Lachrymal Probe
 Cylindrical Tip, Malleable S.Steel
 B#40-230 1/2
 B#40-330 0/00
 B#40-430 000/0000
 scale : 1



Crawford Retrieving Hook
 B# 40-600
 scale : 1

Lachrymal Dilator
 B#40-500
 scale : 1



Crawford Design
 120mm long lachrymal probes attached to
 300mm long silicon tube
 Single use sterilized by ETHYLENE OXIDE



D.C.R. Bodkin Design
 Mini lachrymal Probes attached to silicon tube
 Single use. Sterilized by ETHYLENE OXIDE

