



1



2 Bina Eye Speculum
Fenestrated open blades
B# 60-100



3 Blade Size
Premature 3mm
Infant Size 5mm
Child Size 10mm
Adult Size 15mm



Lieberman Lid Speculum(K-wire)
B# 60-300
Blade Size 14mm



Lieberman Lid Speculum With
Aspirating Cannula
B# 60-350
Blade Size 14mm



Reversible Speculum
Solid Plates
B# 60-340
Plate Size 13.5mm



Eye Speculum
Spring, Solid Plates
B# 60-100-p
Plate Size 13.5mm



Lieberman Speculum
Solid Plates
B# 60-320
Plate Size 13.5mm



Sauer Speculum
Open Plates
B# 60-50
Plate Size 8mm